

IMASONIC reference: cdc9420-8 dated July 1st 2021

Liability

In the context of this proposal and any future studies and achievements, the Customer agrees that it is his responsibility to check, or to have checked, before using the device on the human body or in any way presenting risks for the human body, that the whole system (including the transducer, the electronics, the order processor, the mechanical positioning elements, etc) functions in a way that is safe for the patient and the user. This responsibility includes taking account of the statutory requirements of the countries concerned and the translation of these regulations into requirements applicable to transducers, taking into account the specificities of the system (for example, the electric structure of the system: insulation and grounding choices, etc) and the integration and operating conditions of the transducer.

IMASONIC shall not be liable for:

- The performance of the system, in particular possible failures of the system or the consequences of such failures to patients,
- The system's non-compliance with legal requirements through insufficient information or an incorrect translation by the Customer of the legal requirements into requirements applicable to transducers,
- The consequences of any use of the transducers that does not comply with the using conditions defined in the present document.

In this context, IMASONIC proposes that its responsibilities are as follows:

- To analyze the requests and needs expressed by the customer and to seek solutions for these requirements,
- To make IMASONIC's expertise and resources available to meet requirements and to guarantee compliance with the agreed specifications,
- To ensure the appropriateness and the performance of the methods and manufacturing processes used in executing the orders,
- To alert the customer if a doubt arises concerning the relevance or the reliability of a study or a manufacturing process with a potential link with the safety of the device.

The current proposal takes into account the customer's requirements that IMASONIC is aware of to date.

Other actions may be proposed by IMASONIC, if required and on the customer's request, for additional projects or to follow up future developments.

OVERVIEW - APPLICATION

The current document presents the specifications of a 256-element array transducer that will be incorporated into a device developed by the customer. This device is also composed of an optical source, a mechanical positioning system, a computing system for the image construction and the user interface.

The device will be used for experimental developments and not yet a clinical stage.

NORMATIVE ASPECTS

Considering the experimental laboratory use without any use on the human body, we have considered that the system including the transducer is not required to comply with any standard. Hence, no specific action has been taken into account neither in the design nor in the testing plan.

Normative requirements can be taken into account on request (see also § "Liability").

IMASONIC reference: cdc9420-8 dated July 1st 2021

Technical Characteristics

ACOUSTICAL CHARACTERISTICS:

Array type:	2D array of square elements distributed over a circular aperture with center hole
Mechanical focusing:	spherical, 40 mm (+1/-2 mm) radius of curvature
Element number:	256
Element dimensions:	See table in Appendix 2
Element numbering:	According drawing in Appendix 3
Centre frequency (-6dB):	4 MHz \pm 10% (*) (**)
Acoustical impedance matching:	water (1.5 MRayl)
Bandwidth (-6dB):	\geq 45 % in Transmit/Receive mode (*) (**)
Homogeneity in sensibility:	\pm 5 dB in T/R mode (*) (**) <i>Taking into account only the elements with equal area</i>
Number of non working elements:	6 maximum (adjacent or not) <i>These elements will not be taken into account in the determination of the mean values of the parameters.</i>

(*) average value for all the measured elements

(**) measured in conditions described in "Final control" chapter.

Note: the transducer structure has been chosen to be a trade-off between bandwidth and sensitivity. However, regarding the difficulty to predict the sensitivity, no commitment can be taken on its value. Consequently, the transducer design shall be validated by the customer regarding the expected performance and image.

MECHANICAL CHARACTERISTICS

Material:	Stainless steel
Geometry and dimensions:	According to drawing 9420D0100000 in Appendix
Output cable:	2 outputs, on rear side of transducer
Probe holder interface:	According to drawing 9420D0100000 in Appendix

INTERCONNECTIONS

Cable:	2 x multi-coaxial cables with overall shielding (85 Ω)
Cable protection:	none (external jacket: PVC)
Cable length:	1.5 meters \pm 0.1m
Connector and wiring diagram:	2x ITT DLM5-260 compatible VERASONICS (see appendix for details)
Ground connection:	Shielding of the cables, shell of the connectors, transducer's housing are all connected together and not connected to the electrical ground.

IMASONIC reference: cdc9420-8 dated July 1st 2021

Recommended using conditions****

Integration of the transducer in the system

- Avoid any mechanical shocks on the front face
- **Avoid any mechanical stress on housing and cable bundle when integrating the transducer in the mechanical set-up.**
- **Avoid any mechanical stress on internal tube. This could damage the active part.**
- Do not handle the probe by the cables
- Do not strain the cable more than the radius of curvature (ROC) below
 - Minimum static ROC: 32 mm
 - Minimum Dynamic ROC: 64 mm

Storage, transport and maintenance

- Transducer is fragile, manipulate with care.
- Storage temperature : 0°C to +50°C
- Storage humidity : 10 to 75%
- Cleaning conditions: wiped with a tissue and 75° isopropyl alcohol taking care of not wearing the front face with particles between the front face and the tissue
- The box supplied with the transducer shall be used for any transportation of the transducer.

Utilization

- Continuous immersion in water
- Pressure: under 0.5 m of water
- Temperature : 10 to 40°C
- Excitation: The transducer is used in receive mode mainly and detection sensitivity and smooth band are the parameters to be privileged in the acoustical design of the transducer.
- The transducer will be connected to a Verasonics electronics (input impedance adjustable from 115 Ohms to 3kOhms)

(****) Other using conditions are not taken into account by IMASONIC

Engraving

On housing

4 MHz – 256 elts

Serial number

Index marking (see drawing in Appendix 1 for description)

On connectors

1→ 128 (for connector 1), 129→256 (for connector 2)

Final Control

Temporal response & frequency spectrum measurement (IMASONIC procedure INP413)

Performed in water in T/R mode, on a concentric spherical target

Excitation signal: 1 negative square, T/2, at the nominal frequency

Warranty

12 months in the above-mentioned using conditions.

APPENDIX 1


General Drawing



IMASONIC reference: cdc9420-8 dated July 1st 2021


Appendix 2 Elements dimensions

Dimensions from the front-face point of view					



IMASONIC reference: cdc9420-8 dated July 1st 2021

Appendix 3 Element Numbering



View from transducer rear face

IMASONIC reference: cdc9420-8 dated July 1st 2021

Appendix 4 Pinout diagram

Connector # 1 : elements 1 to 128
Connector #2 : elements 129 to 256

Pin	1	2	3	4	5	6	7	8	9	10	Pin
A	GND	GND						GND	GND	GND	A
B	GND									GND	B
C	GND									GND	C
D	GND									GND	D
E	GND									GND	E
F	GND									GND	F
G	GND									GND	G
H	GND									GND	H
J	GND									GND	J
K	GND	EI 93	EI 94	EI 95	EI 96	EI 97	EI 98	EI 99	EI 100	GND	K
L	GND	EI 92	EI 91	EI 90	EI 89	EI 104	EI 103	EI 102	EI 101	GND	L
M	GND	EI 85	EI 86	EI 87	EI 88	EI 105	EI 106	EI 107	EI 108	GND	M
N	GND	EI 84	EI 83	EI 82	EI 81	EI 112	EI 111	EI 110	EI 109	GND	N
Space for Connector Latch Mechanism											
P	GND	EI 77	EI 78	EI 79	EI 80	EI 113	EI 114	EI 115	EI 116	GND	P
R	GND	EI 76	EI 75	EI 74	EI 73	EI 120	EI 119	EI 118	EI 117	GND	R
S	GND	EI 69	EI 70	EI 71	EI 72	EI 121	EI 122	EI 123	EI 124	GND	S
T	GND	EI 68	EI 67	EI 66	EI 65	EI 128	EI 127	EI 126	EI 125	GND	T
U	GND	EI 4	EI 3	EI 2	EI 1	EI 64	EI 63	EI 62	EI 61	GND	U
V	GND	EI 5	EI 6	EI 7	EI 8	EI 57	EI 58	EI 59	EI 60	GND	V
W	GND	EI 12	EI 11	EI 10	EI 9	EI 56	EI 55	EI 54	EI 53	GND	W
X	GND	EI 13	EI 14	EI 15	EI 16	EI 49	EI 50	EI 51	EI 52	GND	X
Y	GND	EI 20	EI 19	EI 18	EI 17	EI 48	EI 47	EI 46	EI 45	GND	Y
Z	GND	EI 21	EI 22	EI 23	EI 24	EI 41	EI 42	EI 43	EI 44	GND	Z
AA	GND	EI 28	EI 27	EI 26	EI 25	EI 40	EI 39	EI 38	EI 37	GND	AA
BB	GND	EI 29	EI 30	EI 31	EI 32	EI 33	EI 34	EI 35	EI 36	GND	BB
CC	GND	GND						GND	GND	GND	CC
Pin	1	2	3	4	5	6	7	8	9	10	Pin